
IFU-023 In03

NOTICE D’INSTRUCTIONS

IMPLANTS DENTAIRES / VIS DE COUVERTURE IMPLANT / PILIERS / FAUX MOIGNONS / ATTACHEMENTS

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INSTRUCTIONS FOR USE

DENTAL IMPLANTS / IMPLANT COVER SCREWS ABUTMENTS / ABUTMENT POSTS / ATTACHMENTS

Edition December 2024

These medical devices are to be placed only by qualified health professionals trained in their use who have read these instructions. Prior formal training is strongly recommended so as to become familiar with the procedures specific to oral implantology, and in particular with Victory® implant systems.

Scope of application :

These instructions apply solely to the dental implant systems sold by Victory®, the implant cover screws supplied with Victory® implants and the sterile abutments, abutment posts and attachments sold by Victory®.

Product description :

These devices are manufactured of titanium (as per ISO 5832-2 or ISO 5832-3).

Indications :

Products intended to serve as artificial dental roots allowing attachment of a prosthesis. Victory® implants, in conjunction with Victory® prosthetic components, can be used for all types of prosthoses, from single-tooth replacements to full-arch rehabilitations, and for the stabilization of removable prosthoses.

Contraindications :

Numerous factors can compromise the result of implant treatment, and in particular osseointegration of the implant and/or the life of the prosthetic restoration. The following factors must always be evaluated before envisaging implant treatment (non-exhaustive list):
 - allergy or intolerance to the materials (suspected or proven) ;
 - insufficient bone quality and/or volume ;
 - oro-dental pathologies ;
 - poor oral hygiene ;
 - para-functional habits, such as bruxing, or excessive masticatory force ;
 - tobacco, alcohol or substance abuse ;
 - psychological impairment ;
 - infection or inflammation (chronic or acute) ;
 - scheduled or ongoing incompatible treatments such as chemotherapy, radiotherapy or bisphosphonates ;
 - systemic pathologies such as diabetes, metabolic bone diseases or immune system disorders, hematological and/or clotting disorders.
 For young patients, it is preferable to wait until they have completed their growth.

Warnings and precautions :

Prior to any implant treatment, the dentist or doctor must make sure that the patient is apt to undergo such treatment. A thorough preimplant work-up is essential (medical history, clinical examination, radiological studies, study models in occlusion, surgical guide).

Information of the patient :

The health professional must give the patient all necessary information concerning the proposed treatment and obtain his/her consent. This information includes the contraindications, possible adverse events, precautions to be taken, and the obligation to undergo postoperative check-ups to prevent any unwanted evolution such as mobility, screw loosening, or excessive bone loss. The patient must inform his/her dentist or doctor immediately of any unusual event such as pain, mobility or loss of the prosthesis or a retaining screw, infection, suppuration around the site of implantation. The patient must systematically inform the necessary persons during any examination (such as MRI) or ulterior treatment that he/she has one or more dental implants.

Complications and potential adverse events :

- Possible peri-operative complications include, but are not limited to: error in implant selection or positioning, implant placement beyond the predefined limit of the osteotomy that risks damaging a critical anatomic element, insufficient depth of drilling that prevents placement of the implant, over-dimensioning of the osteotomy that prevents adequate retention of the implant, deformation of the implant mount, difficulty for removal of the implant mount screw, placement of a cover screw off axis that is susceptible to damage the implant threads, perforation, fracture, nerve injury.
- Possible transitory postoperative events include, but are not limited to: edema, hematoma, infection, pain, haemorrhage, soft tissue inflammation, altered sensitivity (hyperesthesia, hyposthesia, paresthesia), speech problems, dehiscence, delayed healing, loosening of prosthetic retaining screws.
- Possible adverse events at long term include, but are not limited to: screw loosening, prosthetic loosening, altered sensitivity, bone or implant loss, fistula, damage to an adjacent tooth, neurological lesion, infection, gum retraction, hyperplasia, esthetic problems, oroantral or oronasal communication, perforation of the vestibular bone plate, diapneusia, periapical granuloma.

Sterility :

Victory® implants are sold sterile; a sterile cover screw is included in the packaging. These products are intended for single-use only and must not be reused. Victory® takes no responsibility in case of reuse (risk of infection and/ or cross-contamination).

Storage :

These devices must be stored in a clean, dry environment, in their original packaging, away from direct sunlight. Do not use a product if the expiration date on the label has gone by or if the packaging has been opened or damaged.

Instruments for implant placement :

Use only dedicated Victory® instrumentation for the placement of these devices. Victory® takes no responsibility in case of use with ancillaries of other origins. The Surgery Manuel describes the protocols recommended by Victory®. Prior to use, make sure that the necessary instruments are available, functional and sterile. In particular, verify that all required instruments attach correctly to the handpiece or high-speed handpiece. Utilize appropriate irrigation during osteotomy to avoid thermal injury to the implant site.
 - Drills are single-use only.
 - Titanium cutters are reusable. Before use in the oral cavity, make sure that the cutter rotates properly when attached to the high-speed handpiece or handpiece. The cutting power of the cutter must be sufficient to limit heat build-up during use. The number of utilizations of a given cutter (maximum 10) depends on wear and any damage. After each use, cutters must be cleaned according to the usual protocols applied in dental practice, disinfected, and sterilized in an autoclave (2.1 bars for at least 18 minutes at 134°C / 274°F applying a validated protocol). For more information, refer to the instructions for the retreatment of reusable Victory® dental instruments.

Prosthetic restorations :

Only prosthetic components sold by Victory® must be used for prosthetic restorations on Victory® implants. The Victory® product catalog details the available products. Cleaning and sterilization are required for all prosthetic components and instruments used during the various prosthetic phases: utilize the usual protocols used in dental practice (sterilization in an autoclave at 2.1 bars for at least 18 minutes at 134°C / 274°F using a validated protocol). Any mechanical overload on an implant is a potential source of failure. Occlusal equilibration and adjustment of the prosthesis are essential.

Responsibility of the user :

The user remains responsible for the treatment of his/her patient and the utilization of these products. The health professional must respect the rules of dental surgery in general and pay particular attention to asepsis. Users are advised to consult the instructions for the placement of Victory® implants and the use of Victory® prosthetic components.

In no case and in no manner, neither Victory® nor its consultants can be considered, directly or indirectly, coauthor of a therapeutic decision or act. The user remains entirely responsible vis-à-vis his/her patient concerning a therapeutic decision to use a given product and/or technique. The health professional must know when unfavorable conditions make it unadvisable to place an implant and when the level of difficulty exceeds his/her degree of competency. Victory® takes no responsibility whatsoever, direct or indirect, in case of an error in judgment or improper technique concerning the use of these products, their ancillaries or associated components.

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