IFU-025 In43

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ANCILLAIRES

CE 0459

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To prevent any risk of cross-contamination, reusable dental instruments marketed by Victory® must be cleaned, disinfected and sterilised before the first use and after each use. It is strictly forbidden to reuse devices marked «single use». Cleaning, disinfection and sterilisation procedures must be carried out by trained personnel who are familiar with this instruction leaflet. It is recommended to wear appropriate protective equipment.

### Scope of application:

This leaflet applies to the reusable ancillary equipment marketed by Victory® for use with its implant systems.

### Description of the device:

**Cutting instruments:**

Stainless steel implant drills and gingival cutters for use with a dental contra-angle. The devices comply with the requirements of ASTM F899.

Titanium alloy cutters that meet the requirements of ISO 5832-3

**Fitting instruments:** stainless steel ASTM F899 and/or titanium alloy that meet the requirements of ISO 5832-3

### Indications :

**Cutting instruments:** Drills designed for maxillary and mandibular osteotomies to fit a Victory® axial dental implant to which a dental prosthesis can be attached. Gingival cutters designed to remove gingival tissue prior to osteotomy.

Cutters designed for maxillary and mandibular osteotomies to fit a Victory® basal dental implant to which a dental prosthesis can be attached.

### -Selection of the drill:

After determining the width and depth of available bone, select the Fractal® drill (FE, FR or FD series) according to the bone density (Misch and Judy classification ). The pilot/pointer drill is used to start the osteotomy.

Pilot drill and pointer [FP-5; FP-15; FPM-32; FPIL; FPiXL]: 1200 rpm

Terminal drill [FE; FD; FR]: 800 to 1000 rpm

Cortical drills: 600 rpm

### -Caution:

For safety, leave a 2 mm gap above the inferior alveolar nerve.

- Before drilling, ensure that the bone crest is clear. Use the tip of the pilot drill to mark the point of the bone to be drilled (optional).

- low bone density: preliminary osteogenic preparation with a manual and/or rotary Victory® Osteotensor® is recommended 45 to 60 days before fitting the implant.

**-Irrigation:** Ensure that there is adequate external irrigation during the osteotomy to avoid thermal damage to the bone.

### Fitting instruments:

devices for tightening and/or loosening dental implants or/and prosthetic components for fitting.

Fitting instrument: 30 rpm

Refer to the Victory® Surgery Manuals for instructions for use.

### Limitations and restrictions on repeated treatment:

The lifespan of the instruments is normally determined by wear and damage resulting from their use. Inspect the device visually before and after each use for signs of wear or damage. Do not use a device that is broken, defective, or suspected to be so. Dispose of any damaged medical devices; follow the applicable legislation for the disposal of medical waste.

Bone density	Number of uses
D1	5
D2-D3-D4	10

### Regulations:

The protocol for repeated treatment must ensure the sterility of the treated device. Always comply with the regulations in force. Consult the following documents:

- Good Hospital Pharmacy Practices, Guide to Good Disinfection Practices, Guide to Good Sterilisation Practices, Guide to the Control of Treatments for Reusable Medical Devices (ref FD 598-135 of April 2005).

- List of dental disinfectants published by SFHH and ADF (2018).

### Dismantlable instruments (fitting instruments):

Any instrument with several parts (e.g. mandrel, ring/transfer holder, ratchet wrench) must be dismantled before cleaning, disinfecting and sterilising. Each part of the instrument must be pre-disinfected, cleaned and disinfected separately. The same applies to instruments assembled for the purpose of the procedure (e.g. slotted screwdriver mounted on a ratchet wrench).

### Choice of instruments:

this depends on the purpose, the use and the clinical case. For more information please refer to the surgical manuals.

### Recommendations for cleaning and sterilisation:

If these recommendations are not followed, this may cause a risk to the patient and the device may reach the end of its lifespan more rapidly.

This cleaning and sterilisation protocol can only be performed by staff with proper training and protection, in accordance with the regulations in force. It is imperative to wear appropriate clothing (mask, gloves and goggles), in order to avoid any risks of infection and injury.

### Cleaning:

Remove new non-sterile devices from their packaging before treatment.

- Immerse the devices in a suitable disinfectant solution. Follow the recommendations of the disinfectant manufacturer for dosage, exposure time and temperature.

- Rinse the devices thoroughly with water for at least 1 minute to minimise the chemical residue remaining on the product.

- Dry each device thoroughly with lint-free sterile drapes (disposable towels).

- Place the device in a sealed sterilisation pouch that conforms to ISO 11607 for sterilisation in an autoclave.

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-Victory® recommends use of a class B autoclave according to EN 13060 for all devices bearing the logo

- The instructions and recommendations of the autoclave manufacturer must be followed.

- The autoclave sterilisation cycle is performed at 134°C for 18 minutes.

### Packaging:

After drying and visual inspection, place the instruments in a pouch. Use sterilisation pouches that comply with current standards and regulations and are suitable for the devices concerned. Each part of a dismantlable instrument must be sterilised separately. Follow the instructions for sealing the pouches.

### Storage:

These devices must be stored in a clean, dry place; do not expose to direct sunlight.

- It is the practitioner's responsibility to define the device's period of sterility, taking into account the type of packaging used and the storage conditions.

### Practitioner's responsibility:

The practitioner is responsible for the treatment of their patient and the use of these devices. They must be sure to comply with general dental surgery regulations and particularly ensure compliance with the principles of sterile technique. The practitioner should refer to the Victory® Implant Surgical Technique before using the devices.

Under no circumstances and in no instances shall Victory® or its consultants be considered, directly or indirectly, jointly responsible for a decision or a therapeutic practice. The practitioner remains fully responsible towards the patient for their therapeutic decision to use a particular product and/or technique. The healthcare professional is responsible for refraining from fitting an implant or an osteosynthesis screw if the conditions are unfavourable. Victory® accepts no direct or indirect responsibility for any errors of judgement or practice when using these devices, ancillary equipment or any associated components.

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