VICTORY
ByVisy

FR

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These devices may only be used by qualified persons who have received training on their use, and who have read these instructions. Prior training is highly recommended to become familiar with the specific procedures involved in oral implantology in general, and the Victory® implant range in particular. The listed devices are sold as non-sterile and single use. Victory® accepts no responsibility for the reuse of a previously implanted device.

Scope:
These instructions apply to all non-sterile prosthetic components with or without surface coating, TiN (Titanium nitride) and/or DLC (Diamond-Like Carbon).

Description of device:
Prosthetic components are connected directly to the endosseous dental implant designed to facilitate prosthetic restoration. Specific features of VOA®: Victory® osteosynthesis screws are used in combination with a plate-form Diskimplant® and/or a Diskimplant® to fix or increase their primary stability.

Indications:
Prosthetic components are devices whose purpose is to create with implants (in the maxilla or mandible) stable, resistant, effective, non-irritogenic and long-lasting anchors to which a removable or fixed prosthesis is fitted to provide patients with no teeth or with missing teeth adequate function, satisfactory comfort and a socially acceptable appearance.

For self-drilling osteosynthesis screws:
These devices are designed to fix a plate-form Diskimplant® in and/or on the maxilla or mandible bone and increase the bone blockage of a Diskimplant®, if necessary. The surgeon must check that the Victory® osteosynthesis screw is properly attached to the dedicated instrument used (Victory® screwdriver reference TMCVO). If the bone is very dense, make a preliminary ø 1 mm hole to a depth of 3 mm.

Contraindications:
Many factors can compromise the outcome of implant surgery including the osseointegration of the implant and/or life of the prosthetic restoration. The following factors should always be assessed before considering implant therapy (non-exhaustive list):

- suspected or known allergy or intolerance to the material- inadequate bone quality and/or volume
- oral/dental diseases- poor oral hygiene
- parafunctions (bruxism) or excessive chewing- smoking, alcoholism, drug use- psychological disorders
- chronic or acute infectious or inflammatory condition- scheduled or ongoing incompatible treatment such as chemotherapy, radiation therapy or IV bisphosphonates
- systemic diseases such as diabetes, conditions affecting bone tissue or the immune system, haematological and/or coagulation problems. In all cases, it is best to wait until the patient has finished growing.

Warnings and precautions:
Before performing any implant treatment, the surgeon must check that the patient is suitable for the implant. A full pre-implant assessment must be carried out (history, x-rays, occlusal study models, surgical guide).

Victory® osteosynthesis screws may only be used in combination with plaque-form Diskimplants® and Victory® Diskimplants®. The Victory® Product Catalogue details the entire product range. The surgery manual describes the protocols recommended by Victory®.

Information for patients:
Surgeons must provide the patient with all necessary information on the proposed procedure and collect the consent of the patient. This information includes contraindications, possible side effects, precautions to be taken, the requirement for post-surgical follow-up to prevent any undesired effects (e.g. mobility, loosening of screw or excessive bone loss). The patient should inform the surgeon immediately of any unusual event such as pain, mobility or loss of the prosthesis or fixation screw, or infection or suppuration around the implant site. Before any investigation (such as an MRI) or subsequent treatments, the patient must always inform the technician/doctor that he/she has one or more implants.

Potential complications and adverse effects:
- Possible per-surgical complications include: incorrect selection or positioning of the osteosynthesis screw, inadequate blocking of the screw, off-axis screwing potentially damaging the implant thread, unscrewing, loosening and gingival retraction.
- Possible treatment post-surgical effects include: swelling, bruising, infection, pain, bleeding, soft-tissue inflammation, sensitivity and/or speech problems, dehiscence, delayed healing, unscrewing of the prosthesis.

- Possible long-term side effects include: change in fit (loosening) or migration of the device which may require further surgery for removal or replacement, sensitivity problems, bone loss, fistula, injury to adjacent teeth, nerve damage, infection, gingival retraction, hyperplasia, cosmetic issues.

Prosthetic restorations:
Only prosthetic components sold by Victory® should be used for prosthetic restorations on Victory® implants. The Victory® Product Catalogue details all of the ranges. Always select prosthetic components based on the emergence to be fitted, the type of prosthesis planned (single tooth, bridge, removable) and the actual clinical case. Balancing the prosthesis and adjusting the occlusion are essential steps as any mechanical overloading of the implant is a potential source of failure and/or of potential damage to the devices.

Cleaning and sterilisation recommendations:
Single-use devices delivered in a non-sterile condition must be cleaned and sterilised before use.

Failure to follow these recommendations may result in a risk for the patient and accelerated wear and tear of the device. This cleaning and sterilisation protocol can only be carried out by correctly trained and protected staff, in accordance with the regulations in force. To avoid any risk of infection and injury, appropriate clothing (mask, gloves and goggles) must be worn.

Cleaning
- Remove the new non-sterile devices from their packaging before treatment.
- Immerse the devices in an appropriate disinfectant solution. Respect the recommendations of the disinfectant manufacturer in terms of dosage, exposure time and temperature.

- Rinse the devices thoroughly with water for at least 1 minute to remove as far as possible any chemicals remaining on the product.
- Dry each device thoroughly with lint-free sterile cloths (disposable towel).
- Place the device in a sterilisation pouch suitable for autoclave sterilisation, sealed according to standard ISO 11607.

Sterilisation :
- Victory® recommends class B autoclave sterilisation according to standard EN 13060 for all devices bearing logo [logo]. The instructions and recommendations of the autoclave manufacturer must be followed. - The autoclave sterilisation cycle is 134°C for 18 minutes.

Storage:
These devices should be stored in a clean and dry place. Do not expose to direct sunlight.

Placement instruments:
Only dedicated Victory® placement instruments with these devices. Victory® accepts no responsibility for the use of equipment manufactured by third parties. Before the procedure, check that the instruments used for the placement is available, operational and sterile.

Responsibility of the practitioner:
The practitioner is responsible for the treatment given to his/her patient and for the use of these devices. He/she shall respect the rules governing dental surgery in general and will pay particular attention to respecting the rules of aseptis. The practitioner should refer to the surgical technique required for Victory® implants before using the devices.

Under no circumstances whatsoever shall either Victory® or its consultants be considered, either directly or indirectly, as being jointly involved in any therapeutic decision or action. The practitioner remains wholly liable vis-à-vis the patient for the therapeutic decision to use a particular product and/or technique. The professional must be able to recognise conditions that are unfavourable to the use of an implant and know when the level of difficulty involved exceeds his or her capabilities. Victory® accepts no responsibility, either directly or indirectly, for any error in judgement or practice when using these devices, ancillary equipment or associated components.

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